Simple Medical(Shenzhen)Limited Instruction for Use

Post-Op Elbow Brace E1

Model: 41003 & 41004 Size: F

1a

2a

3a


3b

3c

4a

5a

5b

Product Structure and Performance

The supports are mainly made of aluminum alloy supports, polymer material accessories, straps, linings and accessories. The product is used to limit flexion and extension and control the range of activities after elbow trauma and repair after surgery. The product adopts streamlined lightweight design. Equipped with a support and comfortable non-slip lining, through the support so that the elbow joint with the upper and lower arms have proper protection, special contour design to provide smooth flow to the elbow and improve patient protection.

INDICATIONS

- Post-operative immobilization or ROM control for elbow.
- Conservative treatment of elbow dislocation or luxation
- The distal humerus, radius and proximal ulna are stable or fixed for fracture

APPLICATION INSTRUCTION

NOTICE: This product must be installed and adjusted by a doctor or professional.

STEP 1. Open the product, remove the strap from the strap buckle, and then position the pivot center of the product link on the external epicondyle of the humerus (external elbow)(Fig 1a).

STEP2: Adjust the length of the connecting rod. Push the telescopic adjustment key on the connecting rod support in the direction of "UNLOCK " to unlock the connecting rod, slide the support to the required length, and then push in the direction of" LOCK " to lock the connecting rod. When wearing, make sure that the telescopic adjustment key of the connecting rod is in the "LOCK " position (Fig 2a).

NOTICE: The lower extension (ulnar side of arm) should extend to styloid process (mass inside wrist). The upper extension (humerus side of the arm) should extend to the armpit or the most comfortable position of the armpit.

STEP3: Fasten the strap. Hold the brace by hand to prevent movement. Start close to the elbow, then wrap the strap down around the arm through the strap, then clip the strap with the crocodile's mouth. The order of ties is as follows:

Strap #1: Fasten the lower elbow.

Strap #2: Fasten upper elbow.

Strap #3: Tie wrist.

Strap #4: Tie the upper part of the arm, don't be too tight. too tight may lead to slow blood circulation or uncomfortable wearing (Fig3a-3c).

STEP 4: Adjust the range of activities. Press the flexion or extension angle adjustment button, move to the marked angle scale as needed, and then press the lock button on the chuck to prevent the connecting rod from moving (figure 4a).

STEP 5: First, the longest strap goes from inside to outside through the "T" strap buckle at the top of the arm. (Fig.5a) The straps were then placed around the neck from the inside to the outside around the wrist through the "T" strap buckle with the alligator's mouth clamped to the straps as needed (5b).

Note:

The locking button on the connecting rod support of the reducer is not in the "UNLOCK" direction for a long time.

Post-operative wearing:

It is up to the doctor to set the angle of hinge brace and the angle of flexion and extension as well as the amount of time to wear each day until your attending physician agrees that you can remove the bracing.

***Please purchase and use the product under the guidance of your doctor or professional.**

* IMPORTANT :

- 1) Products are for individual patients and are not recommended for multiple patients.
- 2) Before adjusting the brace, the patient should consult his own doctor.

3) This product shall be purchased and installed on prescription by a medical practitioner and worn.

4) Proper rehabilitation and activity improvement are also indispensable processes in safe treatment planning.

Strap length adjustment

Determine the required length, open the Velcro at the strap end, separate the Velcro at the strap end stick and strap, cut strap to the desired length, and then stick the Velcro at the strap end back to the end of the strap.

LOCKING INSTRUCTION

- 1) Determine the desired locking angle
- 2) Slide the adjustment button toward the center of the hinge to lock the hinge angle.

CARE AND MAINTENANCE

Liners can use neutral detergent and 30 °C water washed by hand, and natural dry. Do not use clothes dryer or other heat source for drying.

Metal parts should not be washed with water.

WARNINGS:

If you feel pain, swelling, change of sensation, or discomfort when using this product, please stop using it and contact your doctor.

Basic application environment

- 1) Ambient temperature: -15 ~ 60 °C
- 2)Relative humidity:≤ 80%

Product materials

- 1)Metal bracket: aluminum alloy/stainless steel
- 2)plastic bracket: polyamide
- 3)Button / buckle: polyformaldehyde
- 4) Rivet: stainless steel
- 5)Strap/pad: nylon / polyurethane / polyester / EVA/ magic paste

Contraindication

- ⚠ Use with caution for unstable fractures;
- ⚠ Use with caution in osteoporosis;
- ⚠ Use with caution for cognitive impairment;
- ⚠ Use with caution for bony ankylosis;

Storage, transportation and expiry date of products

- 1) Term of validity: 5 years.
- 2) Storage of the product: please put it in a dry and clean room before use to ensure indoor ventilation.
- 3) Transportation: the product shall be dustproof, moisture-proof and pollution-proof during storage and transportation.

Waste disposal

Disposal of wrapping paper (bags): please dispose of the packaging of this product according to the local laws.

Symbol and picture description

- Manufacture Date
- Batch Code
- medical device
- CE Symbol
- Protect from sunlight
- Hand wash only
- Do not use when the packing is damaged
- Warnings and Precautions
- Reference instructions

Accessories list

Instructions ----- 1 copy

After-sales warranty

Warranty: all products from the date of purchase as the original date, normal wear and tear of the product in the use of the process is not considered defective.

The following parts from the date of purchase to the warranty period, if there is any material or process defects, please contact your seller directly, Simple medical will completely or partially repair or replace the braces and accessories.

Metal bracket ---- 1 year

Simple Medical (Shenzhen) Limited
Address: 519 Wisdom Valley, No.1010 Bulong Road,
Minzhi Street, Longhua District,
Shenzhen, China.518109
Tel : +86-755-86038563
Web: www.simplemedical.cn
E-mail: info@simplemedicla.cn`

SUNGO Europe B.V.

Olympisch Stadion 24, 1076DE Amsterdam, Netherlands

***Please refer to the Instructions before usage.**

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