
## Simple Medical(Shenzhen)Limited Instruction For Use

### SIMPLE OA Knee Sleeve X1

Model: 11005&11006 Size: S, M, L, XL, XXL

#### PRODUCT STRUCTURE AND PERFORMANCE

The product is mainly made of high elastic fabric body, metal bracket, straps, key and accessories. Lightweight design with support and soft liner. The unique patented hinge technology can effectively reduce the load of knee varus or valgus. By optimizing the angle adjustment, it helps to relieve the pain and pressure points of the knee joint. The high elastic fabric makes the brace and thigh/calf properly protected, and the special contour design improves the patient's ligament protection and comfort. The product is suitable for a variety of uses, so that patients get a good recovery.

#### INDICATIONS:

- Mild to severe unicompartmental knee osteoarthritis
- Meniscus procedures/repairs
- Medial or lateral compartment

#### APPLICATION INSTRUCTIONS:

**NOTICE:** Before use, please set the hinge of 0 ° on horizontal position. The setting of the support angle should be performed by a trained professional in accordance with the doctor's instructions for inner and outer unloading, which will be adjusted in step 5.

**STEP 1.** Untie two straps and open the brace. Knee bends 30 °, align the position of the hinge center and the center of the knee alignment, and then hold in place to prevent slide (figure 1~3).

**STEP 2.** Pull the four pieces over the knees and legs in turn 1-2-3-4 and fasten them firmly (Fig .4~5).

**STEP 3.** Pull the calf strap #5 through the buckle and tie it tight. Repeat thigh strap #6 (fig.6).

**STEP 4.** Product initial installation: knee bend30 °, turn hinge center screw 1 / 8 (important) clockwise with adjustment tool to increase unloading on the outer side or unclockwise 1 / 8 turn to the center. This is the best starting point for initial use (fig 7 ~8).

**STEP 4-1.**Stand and walk short distance to ensure comfort after each adjustment and support. If there is any discomfort in the knee or leg, adjust the angle slightly back.

If more unloading is required, turn hinge center screws no more than 1/8 turn. Walk after adjustment and check for comfort.

#### \* IMPORTANT:

**Excessive adjustment of hinge center screws may cause discomfort or skin pressure problems. Please adjust carefully**

#### "WEARING TIME" :

Gradually increase the wearing time of the Brace and adjust the wearing time of the brace if your body allows.

#### PREOPERATIVE WEARING TIME:

The first wearing time starts from 30 minutes and gradually increases to no less than 2 hours per day according to the degree of body adaptation. The supporting angle of hinge can be adjusted again in condition of no discomfort.

#### POSTOPERATIVE WEARING:

The standard of wearing should be decided by the doctor, including the angle of hinge support, the angle of flexion and extension, and the time of wearing every day.


The pain caused by osteoarthritis takes time to subside. Wearing a brace regularly for weight-bearing activities (such as walking) will start to feel lighter as it will reduce joint load. This may require wearing it for weeks or even more time to reduce joint load until your primary care physician allows you to remove the brace.

#### RANGE OF MOTION ADJUSTMENTS

Start from the standard 0 °, additional flexion and extension angle accessories can be installed and they shall be fitted by the doctor as required.

Optional extension angles include 5°, 10°, 15°, 20°, and 25°; Optional flexion angles include 30°, 45°, 60°, 75°, and 90°.

\*Please purchase and use the product under the guidance of your doctor or professional.

#### CARE AND MAINTENANCE

High elastic fabric body can use neutral detergent and 30 °C water washing by hand, and natural dry. Do not use clothe dryer or other heat source for drying.

Metal parts should not be washed with water. Hinges can be lubricated with dry silicone spray.

#### WARNINGS:

If you experience pain, swelling, sensory changes, or discomfort with this product, please stop using and contact your doctor.

#### FOR SINGLE PATIENT USE ONLY

**NOTICE:** The device is a one-element holistic treatment program managed by medical professionals. There is no guarantee that damage can be prevented by using the device.

#### BASIC PRODUCT ENVIRONMENT

- 1) Ambient temperature: -15 ~ 60°C
- 2) Relative humidity: ≤80%

#### PRODUCT MATERIALS

- 1) High elastic fabric body: spandex fiber
- 2) Metal support: aluminum alloy/stainless steel
- 3) Buckle: nylon
- 4) Angle adjusting key: nylon/carbon steel
- 5) Strap: nylon/polyurethane/polyester

#### CONTRAINDICATION

- ⚠ Use with caution for unstable fractures;
- ⚠ Use with caution in osteoporosis;
- ⚠ Use with caution for cognitive impairment;
- ⚠ Use with caution for bony ankylosis;

#### STORAGE, TRANSPORTATION AND EXPIRY DATE OF PRODUCTS

- 1) Term of validity: 5 years.
- 2) Storage of the product: please put it in a dry and clean room before use to ensure indoor ventilation.

3) Transportation: the product shall be dustproof, Moisture-proof and pollution-proof during storage and transportation.

#### WASTE DISPOSAL

Disposal of wrapping paper (bags): Please dispose of the packaging of this product according to the local laws.

#### SYMBOL AND PICTURE DESCRIPTION

Manufacture Date

Batch Code

medical device

It is valid for five years

Protect from sunlight

Hand wash only

Warnings and Precautions

Reference instructions

#### ACCESSORIES LIST

- Instructions ----- 1 COPY
- Angle adjusting accessories ---- 1 SET
- Angle adjustment key ---- 1PCS
- Hexagon wrench ---- -- 1PCS
- Knee pad ---- - 1 PCS (sets)

#### AFTER-SALES WARRANTY

Warranty: all products from the date of purchase as the original date, normal wear and tear of the product in the use of the process is not considered defective.

The following parts from the date of purchase to the warranty period, if there is any material or process defects, please contact your seller directly, simple medical will complete or partial repair or replacement of equipment and accessories.

Metal support ---- -- 1 year

**Simple Medical (Shenzhen) Limited**  
Address: 519 Wisdom Valley, No.1010 Bulong Road,

Minzhi Street,Longhua District,  
Shenzhen, China.518131

Tel : +86-755-86038563

Web: www.simplemedical.cn

E-mail: service@simplegroup.net

SUNGO Europe B.V.

Fascinatio Boulevard 522 ,Unit1.7,2909VA Capelle aan den IJssel ,The Netherlands

Tel : +31(0)10-3034500

\*Please refer to the Instructions before usage

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